

## Appendix: Mechanical and Leakage Integrity Testing Protocols for Evaluating the Performance of Tissue Containment Systems Used During Power Morcellation Procedures

This Appendix lists protocols for the test methods listed in Table 1 to evaluate mechanical strength and leakage potential of the TCS materials.

### Tensile Testing

The tensile testing protocol is based on the following standards:

Standard	Aspects of Standard Used/Incorporated
ASTM D412	Section 10.1: Die – Type C Dye
ASTM D882	All with the Exception of Section 6

The majority of the test procedure is derived from ASTM D882. However, in place of section 6 of ASTM D882 this standard, test samples should be cut using an ASTM D412 Type C (25 mm × 115 mm) dumbbell shaped die.

Additional information for the tensile testing can be obtained from [Herman et al. \[2023\]](#).

### Burst Testing

Information for the burst testing can be obtained from [Herman et al. \[2023\]](#).

SolidWorks® files of the test rig can be downloaded [here](#).


### Puncture Testing

Resistance to puncture for each containment system was obtained by measuring the force required to cause a standardized puncture pin to completely penetrate through the thickness of the specimen. A 50mm circular coupon is sandwiched between two circular specimen holders with an internal test diameter of 25 mm using four screws (Figure 1). The specimen holder is then secured to a puncture test fixture (Figure 2) which is




designed to be secured on top of a 50 mm compression plate. For the puncture pins, one of the two standard durometers pins, Type OO and Type D from ASTM D2240 are used. The pin is moved uniformly downwards with a speed of 25 mm/min until it punctured through the thickness of the containment system while the force-displacement data was recorded. The threshold force at which the pin traversed through the full thickness of the TCS (i.e., fully punctured) was called the full-puncture force ( $F_{\text{puncture}}$ ).


Additional Information for the puncture testing can be obtained from [Herman et al. \[2023\]](#). SolidWorks® files of the sample sandwich, puncture rig stand, and bracket can be found [here](#).



\* Same for all Configurations



Top Cap




Bottom Cap

Units Are mm unless Specified or Dual Dimensioned

Aluminum Caps (Top and Bottom) and 2 Rubber Gaskets for 1 assembly  
- Caps can be made from 2" OD, 1" ID Round Stock


Gasket should not be thicker than 3/32" (Inch)

Can use counterbore instead of countersunk as long as screws provided are flush with both sides





**Figure 1: CAD Drawing of the Sandwich Assembly Used in the Puncture Test and Partial Puncture/Dye Test Setup**



**Figure 2: CAD Drawing of the Puncture Rig Setup Base**

## Dye Testing

The protocol for dye testing has been adopted from ASTM F1670. All sections with the exception of Section 7: Reagents from ASTM F1670 are applicable here. Annex A of ISO 16603 may be used in place of Section 7 from ASTM F1670 for artificial blood soil. [Herman et al. \[2020\]](#) also provides additional information about dye testing.




## Bacteriophage Testing


The protocol for bacteriophage testing has been adopted from the protocol outlined in ASTM F1671. All sections are applicable for the Bacteriophage Testing procedure for this RST. [Herman et al. \[2020\]](#) also provides additional information about bacteriophage testing.

## Partial Puncture / Dye Testing

Information for the partial puncture testing can be obtained from [Herman et al. \[2023\]](#). SolidWorks® files of the modified dye testing cell, and mounting plate can be found [here](#).



**Figure 3: CAD Drawing of the Modified Dye Testing Cell to be Used in the Partial Puncture and Dye Testing.**



**Figure 4: CAD Drawing of the Mounting Bracket for Modified Dye Testing Cell**